

Ovens & Furnaces



Despatch Double Oven: up to 550°F (consistent)

We could accept cement cores for drying, but we'd have to open a hole in the ceiling to drop them in. That is what this oven was used for by the BOR, before we gave it a home here in Worland. We use this oven for drying multiple samples, large samples, whopping big samples..... etc.



Thermolyne Type 1500 Furnace (for melting stuff): maximum 2000+°F

I'm sure I'll understand additional potential uses for this furnace far better after taking the Assay Course at Montana Tech. This is a muffle furnace. It attains a sustained 2000+F, and we plan to use it to heat some ochre samples this summer (2010) at a sustained 650-degrees C.



Thelco Model 18: up to 450°F (consistent)
We use this oven just about every day and the standard setting is at 180-degrees.
We have a microwave. We mostly use it for snacks.....